



#### Incubator with natural air circulation

- Bench top model
- Natural convection
- Entirely made of stainless steel
- Glass door (option Inox door insulated)
- Electronic control
- Door lockup with key
- Heat insulated with fiberglass
- Temperature range from **30°C** to **70°C**
- Time regulation from 0 to 600 minutes and infinite time
- Digital display for temperature and time
- Resolution 0,1°C
- Visual and ringing alarms for :
  - over temperature
  - temperature fall
  - probe damage
  - broken heating element
  - black out
- Voltage 230V AC (110V on demand)
- With solid Schuko plug

#### Incubateur sec avec circulation naturelle d'air

- Modèle de table
- Circulation naturelle
- Entièrement construits en acier inox
- Porte en verre (optionnelle porte inox isolé)
- Contrôle électronique
- Fermeture à clé de la porte
- Isolé par fibre de verre
- Thermorégulateur de **30°C** à **70°C**
- Temps sélectionnable de 0 à 600 minutes et temps infini
- Affichage digital de la température et du temps
- Résolution 0,1°C
- Alarmes visuelles et sonores pour :
  - sur température
  - chute de température
  - échec sonde
  - résistance en panne
  - black out
- Alimentation 230V (110V sur demande)
- Schuko câblé

Art.	Internal dimensions (mm) Dimensions intérieur (mm)	External dimensions (mm) Dimensions extérieur (mm)	Capacity (Liters) Capacité (Litres)	Max power (W) Puissance max (W)	Removable shelves Plateaux extractibles
<b>2431/V</b>	385 x 235 x 218 (h)	560 x 370 x 380 (h)	20	100	2